
SPINE

VERTEBRIS zervikal

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V E T E R I N A R Y

I N D U S T R I A L

VERTEBRIS cervical
Full-endoscopic Spinal
Instrumentation

VERTEBRIS cervical, full-endoscopic techniques

Table of contents

VERTEBRIS cervical	4
Introduction	4
The full-endoscopic posterior technique	6
• Positioning of patient	6
• Determination of access route	6
• Creation of access	7
• Performance of operation	7-9
The full-endoscopic anterior technique	10
• Positioning of patient	10
• Determination of access route	10
• Creation of access	11
• Performance of operation	12-13
VERTEBRIS Instruments	14
• VERTEBRIS cervical posterior basic set	14-15
• VERTEBRIS cervical anterior basic set	16-18
• Accessories	19

VERTEBRIS cervical

Introduction

Posterior access for the full-endoscopic cervical operation

In the area of the cervical spine, radicular symptoms such as pain in the arm are caused by degeneration typically from mediolateral to lateral herniated discs or stenoses of the intervertebral foramen. Such clinical symptoms were first classified as being topographically related to changes in the cervical discs at the beginning of the 1940's. Despite the often good results achieved with conservative procedures, a surgical approach can become necessary in cases of persistent pain or neurological deficits.

The posterior surgical access to the cervical spine was also developed in the 1940's. The anterior access was described at the end of the 1950's. In the meantime, anterior decompression and fusion has become a standard procedure in the surgical treatment of cervical radiculopathy. It is considered a safe and adequate procedure with good fusion rates. Nevertheless, specific problems, such as sintering of the implants, pseudarthroses or access complications

are reported. Adjacent disc degeneration has been discussed as one particular disadvantage of fusion. Here, an attempt is made to reconstruct the intervertebral space with disc prostheses while retaining mobility of the segment. With lateral pathology, the most common alternative to the anterior procedure is the posterior foraminotomy. This is performed without additional stabilisation and therefore allows the segment to retain its mobility. The problems that can occur here include neck pain resulting from the access or intraoperative bleeding. A reconstruction of the intervertebral space is not possible.

In cervical disc herniations with radicular symptoms, the volume of the escaped disc material is usually small. This means that due to the access in the anterior and posterior standard technique, the operation is comparatively extensive relative to the limited pathology. To reduce the disadvantages of the conventional procedure, modifications have been described such as anterior decompression without fusion, various techniques of anterior foraminotomy or posterior microscopically assisted or endoscopically assisted "keyhole foraminotomy". With anterior techniques without reconstruction of the intervertebral space, in particular, possible problems of sintering and segmental kyphosis have been discussed.

Full-endoscopic operations on the cervical spine have been described since the 1990's mainly using an anterior, transdiscal approach. The narrow anatomical situation represented a problem and allowed only the use of small endoscopes and working sleeves. This, in turn, presented technical problems, such as the poor view, working under X-ray control without direct visualisation

or restricted bone resection. Foraminal herniations could not be operated on ventrally.

Anterior access for the full-endoscopic cervical operation

Soft disc herniations are the main indication

With the development of new endoscopes, instrument sets and operating techniques, the technical problems have been overcome and today it is possible to treat cervical herniated discs full-endoscopically and with continual visualisation using an anterior or posterior approach. The possibility of adequate bone resection under vision (for example in the area of the foramen, the uncinat process or posterior edge of the body of a vertebra) and various

operating instruments create technical conditions similar to those of conventional microscopically assisted surgery however with the advantages of a full-endoscopic procedure with 25° telescopes and continuous irrigation.*


The main indications for cervical full-endoscopic operations are "soft" hernia-

Cervical myelon with spinal nerve and disc herniation

ted discs with radicular symptoms, in other words, pain in the arm. Since the cervical myelon must not be manipulated in a medial direction, the posterior access is used in herniations in which the main part is located lateral to the lateral edge of the myelon. Using an anterior approach in this situation, there is no guarantee of reaching the pathology even with resection of the dorsal unciniate process. Herniations located mainly medial to the lateral edge of the myelon count as indications for the anterior access, since here the myelon prevents the posterior approach. For an anterior access, the height of the ventral edge of the intervertebral space in reclinaton must be at least 4 mm to prevent injury resulting from the approach. With either access, possible craniocaudal sequestration must not exceed half the body of the vertebra.

The surgeon must also be capable of the conventional and maximally invasive procedures on the cervical spine. Possible problems and complications of cervical interventions may, in some cases, have serious repercussions. Injury to blood vessels can never be fully excluded in cervical spine surgery and, if this occurs, immediate conversion

The lateral margin of the myelon is the line indicating posterior or anterior access

to an open technique is necessary. Personnel and equipment must always be available for this emergency option.

* Ruetten S, Komp M, Merk H, Godolias G: Full-endoscopic cervical posterior foraminotomy for the operation of lateral disc herniations using 5.9-mm endoscopes: A prospective, randomized, controlled study. *Spine* 2008; 30:940-948

Ruetten S, Komp M, Merk H, Godolias G: A new full-endoscopic technique for cervical posterior foraminotomy in the treatment of lateral disc herniations using 6.9-mm endoscopes: prospective 2-year results of 87 patients. *Minim Invas Neurosurg* 2007;50:219-226

Department of Spine Surgery and Pain Therapy
Head: Dr. Sebastian Ruetten, M.D.

Center for Orthopaedics and Traumatology
St. Anna Hospital, Herne, Germany
Director: Prof. Georgios Godolias, M.D.


VERTEBRIS cervical

The full-endoscopic posterior technique

1. Positioning of patient

The patient is prone with a pelvic and thoracic cushion. For a posterior approach to the cervical spine, lordosis of the head and cervical spine must be straightened and fixed permitting intra-operative X-ray control in two planes. Stabilisation with Mayfield tongs or similar is helpful and allows conversion to an open intervention in an emergency. For the lower cervical spine, in particular, it may be necessary to tape the shoulders in a caudal direction or to apply caudal traction to the arms to relocate the shoulders. A C-arm is required during the operation.

Prone position, head stabilised in Mayfield tongs, caudal traction applied to arms

2. Determination of access route

Using fluoroscopic guidance, the access is selected in the orthograde lateral and posterior anterior view based on the anatomical landmarks. It must be located precisely over the zygapophyseal joints at the level of the disc.

Drawing the line of the zygapophyseal joints in posterior anterior view


Skin incision

Identifying the level of the disc in orthograde lateral view with a cannula and selecting the point of entry

3. Creation of access

After selecting the point of entry through the skin and the incision, the dilator is introduced under lateral fluoroscopic guidance control until it contacts the bone of the zygapophyseal joints. The working sleeve with an oblique opening is inserted over the dilator and the dilator is removed.

Introducing the dilator as far as the zygapophyseal joints

The operating sleeve is introduced over the dilator

4. Performance of operation

The endoscope is introduced through the working sleeve. The operation is performed under vision using different instrument sets via the intraendoscopic working channel and under continuous irrigation. In almost all cases, foraminotomy requires bone resection with various instruments. After dissecting the bony structures, this begins at the descending part of the joint and the cranial lamina depending on the anatomy and pathology. Following this, parts of the

caudal lamina and ascending facet are resected. Even at this stage, care must be taken to avoid spinal nerves and the vertebral artery. The ligamentum flavum is opened and the spinal channel can be entered to resect the herniated disc.

The surgeon works with the endoscope through the operating sleeve

VERTEBRIS cervical

The full-endoscopic posterior technique

Bony parts of the joint and the laminae are resected to open the foramen





Fluoroscopic guidance can be useful while using the burrs and working in the spinal canal

Opened foramen with a view of the ligamentum flavum

View into the spinal channel with the cervical myelon and spinal nerve

After removal of the lateral ligamentum flavum and dissection of the neural structures, the disc herniation can be removed

The sealing caps for the telescope and working sleeve should only be used briefly if bleeding impairs the view. If an operation takes a long time and the blockage of the outflow of irrigation fluid goes unnoticed, effects of volume overload and increased pressure within the spinal channel and neighbouring structures cannot be fully excluded.


Manipulation of the cervical myelon must be avoided at all costs. In general, as with all new techniques, there is an increased risk of complications during the learning phase. This can be further increased due the general anatomical conditions in the area of the cervical spine as compared with the lumbar spine.

VERTEBRIS cervical

The full-endoscopic anterior technique

1. Positioning of patient

The patient is in the supine position. For an anterior approach to the cervical spine, the head and cervical spine must be slightly reclined and stabilised permitting intraoperative X-ray control in two planes. Stabilisation with Mayfield tongs or similar is helpful and allows conversion to an open intervention in an emergency. For the lower cervical spine, in particular, it may be necessary to tape the shoulders in a caudal direction or to apply caudal traction to the arms to relocate the shoulders. A C-arm is required during the operation.


Dorsal position, head stabilised in Mayfield tongs, caudal traction applied to arms

2. Determination of access route

The access is established contralateral to the pathology. The ventral spine is palpated with the fingers while the oesophagus and parts of the trachea are displaced medially and the vessels and nerves laterally. Under fluoroscopic guidance and taking into account the anatomy and pathology, the access is selected precisely over the intervertebral space in the orthograde lateral view.


The ventral spine is palpated

The access is selected over the intervertebral space

3. Creation of access

After selecting the point of entry through the skin and the incision, the first thin dilator is introduced into the intervertebral space under lateral fluoroscopic guidance. Care must be taken to puncture the disc anteriorly and avoid missing it laterally. This prevents not only the subsequent operation but can also lead to injury to the vertebral artery, spinal nerve or oesophagus. As an alternative, the disc can be punctured with a spinal cannula introduced over a guide wire. The first dilator can then be introduced over this. After puncture of the disc with the dilator or spinal cannula, the position is checked in a posterior anterior C-arm view. The remainder of the intervention is performed in the lateral view. The combined dilator sleeve system is introduced into the intervertebral space over the first dilator. The dilators are removed, the operating sleeve remains in the intervertebral space.

Introduction of the thin dilator into the intervertebral space

Introduction of the combined dilator-sleeve system

The operating sleeve remains in the intervertebral space


VERTEBRIS cervical

The full-endoscopic anterior technique

4. Performance of operation

The endoscope is introduced through the working sleeve. The operation is performed under vision using different instrument sets via the intraendoscopic working channel and under continuous irrigation.

To achieve topographic orientation, dissection of the uncinete process, dorsal edge of the body of the vertebra and dorsal anulus is performed. In many cases, bone resection with various instruments is necessary to reach the epidural space. Depending on the anatomy and pathology, the dorsal anulus and posterior longitudinal ligament are opened to gain access to the spinal canal and resect the disc herniation.


Bone resection is often necessary to reach the spinal canal



Working with the endoscope through the operating sleeve

The fluoroscopic guidance can be helpful for orientation during the intervention

Depending on the findings, the posterior longitudinal ligament is opened

Removal of the disc herniation

The sealing caps for the telescope and working sleeve should only be used briefly if bleeding impairs the view. If an operation takes a long time and the blockage of the outflow of irrigation fluid goes unnoticed, effects of volume overload and increased pressure within the spinal canal and neighbouring structures cannot be fully excluded.

Manipulation of the cervical myelon must be avoided at all costs. In general, as with all new techniques, there is an increased risk of complications during the learning phase. This can be further increased due the general anatomical conditions in the area of the cervical spine as compared with the lumbar spine.

VERTEBRIS




Basic set, **VERTEBRIS** cervical posterior
devised by Dr. Ruetten

Article		Types
	Discoscope, 25°, WC 3.1 mm	892108253
	Light cable	8061.353
	Dilator, 2-channel Ø 5.9 mm	8792.764
	Working sleeve, WL 80 mm, 30° bevel	89220.7007
	Working sleeve attachment, Ø 7 mm	89200.1007
	Elevator, Ø 2.5 mm	89250.2025
	Dissector, Ø 2.5 mm	8792.591
	Microrongeur, Ø 2.5 mm, WL 290 mm	89240.2025
	Rongeur, Ø 3 mm, WL 290 mm	89240.3003
	Micropunch, Ø 2.5 mm, WL 290 mm	89240.2225

VERTEBRIS

Basic set, **VERTEBRIS** cervical posterior
devised by Dr. Ruetten

Article	Types
 Bone punch, Ø 2.5 mm, WL 290 mm	89240.2325
 Tube shaft punch, Ø 3 mm, WL 290 mm	89240.3903
 Trigger flex handpiece, complete	8792.6911
 Trigger flex bipolar electrodes (pack of 6)	4792.6912
 Probe with flexible tip, Ø 2.5 mm, WL 290 mm, consisting of:	892506625
 Inner sheath	892506625
 Handpiece	892500600
 Tube sheath	15570644
 Oval burr with lateral guard, Ø 2.5 mm	8792.312
 Oval burr with lateral guard, Ø 3 mm, WL 350 mm	89970.1503
 Oval burr with frontal guard, Ø 3 mm, WL 350 mm	89970.1513



VERTEBRIS

Basic set, **VERTEBRIS** cervical anterior
devised by Dr. Ruetten

Article		Types
	Cervical discoscope, 25°, WL approx. 150 mm	892106250
	Telescope connecting piece	892006000
	Light cable	8061.353
	Working sleeve, 3.8 x 6.2 mm, WL 102 mm	892206038
	Working sleeve, 4.1 x 6.7 mm, WL 102 mm	892206041
	Dilator conical, for working sleeve 892206038	892206438
	Dilator conical, for working sleeve 892206041	892206441
	Dilator, for working sleeve 892206038	892206538
	Dilator, for working sleeve 892206041	892206541
	Handpiece attachment for working sleeve 892206038	892006038
	Handpiece attachment for working sleeve 892206041	892006041

VERTEBRIS

Basic set, **VERTEBRIS** cervical anterior
devised by Dr. Ruetten

Article		Types
	T-handpiece, Ø 12.0 mm	892006120
	Guide rod, Ø 1.8 mm, WL 250 mm	892206318
	Guide cannula, Ø 1.8 mm, IØ 0.9 mm, WL 250 mm	892206118
	Dilatation set, three-part	892206500
	Spinal cannula set, Ø 1.25 mm, WL 90 mm	492206112
	Trepine, Ø 3.6 mm, WL 100 mm, with tissue guard	892606036
	Trepine, Ø 4.0 mm, WL 100 mm, with tissue guard	892606004
	Probe with flexible tip, Ø 2.5 mm, WL 290 mm, consisting of:	892506925
	Inner sheath	892506625
	Handpiece	892500600
	Tube sheath	15570644
	Hooked probe, Ø 2.0 mm, WL 290 mm	892506003

VERTEBRIS

Basic set, **VERTEBRIS** cervical anterior
devised by Dr. Ruetten

Article	Types
 Mikro-rongeur, Ø 2.5 mm, WL 290 mm	89240.2025
 Rongeur, Ø 3 mm, WL 290 mm	89240.3003
 Micropunch, Ø 2.5 mm, WL 290 mm	89240.2225
 Micropunch, Ø 3 mm, WL 290 mm	89240.3023
 Bone punch, Ø 2.5 mm, WL 290 mm	89240.2325
 Bone punch, Ø 3 mm, WL 290 mm	89240.3903
 End-cut burr, Ø 3 mm, WL 350 mm	89260.1113
 Trigger flex handpiece, complete, Ø 2.5 mm	8792.691
 Trigger flex bipolar electrodes (pack of 6)	4792.6912
 Oval burr, Ø 2.5 mm, WL 350 mm	8792.312
 Oval burr with lateral guard, Ø 3 mm, WL 350 mm	89970.1503
 Oval burr with lateral and frontal guard, Ø 3 mm, WL 350 mm	89970.1513

VERTEBRIS

Accessories

Article		Types
	PowerDrive Art1 shaver system, incl. network-cable and CAN-bus connecting cable, 230V, 50/60 Hz	2304.0011
	PowerDrive Art1 shaver system, incl. network-cable and CAN-bus connecting cable, 100V, 50/60 Hz	2304.0021
	PowerDrive Art1 shaver system, incl. network-cable and CAN-bus connecting cable, 110V, 50/60 Hz	2304.0041
	PowerDrive Art1 shaver system, incl. network-cable and CAN-bus connecting cable, 115V, 50/60 Hz	2304.0061
	PowerDrive Art1 shaver system, incl. network-cable and CAN-bus connecting cable, 120V (USA), 50/60 Hz	2304.0071
	PowerDrive Art1 shaver system, incl. network-cable and CAN-bus connecting cable, 127V, 50/60 Hz	2304.00121
	PowerDrive Art1 shaver system, incl. network-cable and CAN-bus connecting cable, 240V, 50/60 Hz	2304.00141
	Power Stick M4 / motor handpiece incl. connecting cable	8564.121
	Double pedal foot-switch	2304.901
	Surgitron radiofrequency unit, 4 Mhz	2343.001/ .002

M E D I C A L

V E T E R I N A R Y

I N D U S T R I A L

VERTEBRIS Lumbar-Thoracic
Full-endoscopic Spinal Instrumentation

VERTEBRIS lumbar-thoracic, full-endoscopic techniques

Table of contents

VERTEBRIS lumbar	4
Introduction	4
The full-endoscopic trans- and extraforaminal technique	6
• Positioning of patient	8
• Determination of lateral access route	8
• Creation of lateral access	9
• Performance of operation	12
• Creation of posterolateral access	13
• Creation of extraforaminal access	14
• Bone resection	15
• Biportal access	16
The full-endoscopic interlaminar technique	17
• Positioning of patient	19
• Determination of access route	19
• Creation of access	20
• Performance of operation	21
• Bone resection	24
VERTEBRIS thoracic	25
Introduction	25
The full-endoscopic transforaminal technique	26
The full-endoscopic interlaminar technique	27
VERTEBRIS basic sets	
• VERTEBRIS lumbar trans-/extraforaminal basic set	28
• VERTEBRIS lumbar and thoracic interlaminar basic set	29
• VERTEBRIS thoracic transforaminal basic set	30
• VERTEBRIS instruments	31-39
Literature	40
Notes	41-42

VERTEBRIS lumbar

Introduction

Lateral approach for the full-endoscopic transforaminal operation

Musculoskeletal pain is one of the most common reasons for visiting the doctor. Degenerative diseases of the spine form part of daily medical practice and their treatment is complicated by medical and socioeconomic problems.

Where severe pain or neurologic deficits persist and all conservative treatment options have been exhausted, surgery may be required. Though conventional operations can achieve good results, damage may ensue due to traumatization. It is therefore important that the surgical techniques used be optimized on a continuous basis.

The latest research results and technical innovations must be critically reviewed so that optimal treatment strategies can be formulated. The aim should be to minimize surgically induced trauma and negative long-term sequelae, taking existing quality standards into account.

Minimally invasive techniques can reduce tissue damage and its consequences. Endoscopic operations with continuous irrigation have advantages that have made them the gold standard in a number of areas. Transforaminal procedures with posterolateral access have been performed in the lumbar spine for more than 20 years now, mostly for intradiscal and intra- and extraforaminal procedures. In our Department of Spine Surgery and Pain Therapy we have therefore been developing a lateral transforaminal and an interlaminar approach for full-endoscopic access to the spinal canal since 1998. These approaches broaden the range of indications and permit the use in selected indications of a visually controlled procedure that is as effective as conventional surgery while benefiting from all the advantages of truly minimally invasive surgery.

Until fairly recently, the endoscopic approach was subject to technical problems in that the intraendoscopic working channel of the available optical systems was small and the repertoire of instruments that could be used was correspondingly limited. Insurmountable difficulties could arise in the resection of hard tissue and in terms of limited anatomic access and mobility. Work on the affected tissue was limited and sometimes had to be performed without direct visualization. New rod lenses with a 4.1 mm intraendoscopic working channel and corresponding new instruments, as well as shavers and burrs, were therefore developed so that full-endoscopic operations could be performed under continuous and precise visual control. This also permitted adequate bone resection. The range of indications for endoscopic spinal surgery was thus

broadened to include intervertebral disc herniations, spinal canal stenosis, and stabilization techniques.

Continuous irrigation provides outstanding intraoperative visualization

The latest generation of endoscopes have a large (4.1 mm) working channel

Full-endoscopic surgery has now won a firm place in the surgery of lumbar spine conditions. Provided that the indications for its use are observed, it represents a useful and safe addition or alternative to conventional surgery. Full-endoscopic operations can also be performed on the cervical and thoracic spine. Recent technical developments and the use of new access routes have led

ling the end of existing operative standards; rather, it should be seen as a valuable additional option within the field of spinal surgery.

Herne, July 2007

Dr. Sebastian Ruetten, M.D.
Head Department of Spine Surgery and Pain Therapy

The development of new instruments broadens the range of possible procedures

to a change that suggests the onset of a revolution in spinal surgery similar to that which occurred in orthopedics with the introduction of arthroscopic procedures. Nevertheless, conventional and maximally invasive operations will continue to play an indispensable role in spinal surgery. Surgeons will need to be able to perform these techniques too in order to overcome problems and complications of full-endoscopic operations such as can occur with any invasive procedure.

The development of full-endoscopic techniques should not be seen as spel-

Department of Spine Surgery and Pain Therapy

Head: Dr. Sebastian Ruetten, M.D.

Center for Orthopaedics and Traumatology

St. Anna Hospital, Herne, Germany

Director: Prof. Georgios Godolias, M.D.

at the Department of Radiology and Microtherapy
University of Witten/Herdecke

VERTEBRIS lumbar

The full-endoscopic trans- and extraforaminal technique

Percutaneous operations aimed at achieving intradiscal decompression of lumbar intervertebral discs were first described in the early 1970s. Optical systems designed for inspecting the intervertebral space after open operation were introduced in the early 1980s. Later, a full-endoscopic technique using a transforaminal approach was developed. In anatomic terms, this means that the intervertebral disc is reached via a posterolateral approach through the intervertebral foramen between the exiting and traversing nerve roots without need for resection of bony or ligamentous structures. The skin entry point for operative access is determined in centimeters from the midline. Most such operations are performed for the purpose of intradiscal or extradiscal foraminal therapy. Reduction of intradiscal volume and pressure can reduce disc-related compression. Removal of intra and extraforaminal disc material is technically possible. Sequestered material lying within the spinal canal can generally be resected in retrograde fashion intradiscally via the annular defect. This is done using an "in-out" technique.

Sequestered nuclear material is found within the spinal canal dorsal to the annulus in the ventral epidural space medial to the medial pedicular line. In many cases it extends to the midline or even to the contralateral side. Clinical experience has shown that the annular defect is often smaller than the diameter of the sequestered material. In addition, there is generally no longer any direct connection to the intradiscal space. In the case of badly degenerated discs or older disc herniations, the continuity of the sequestered material has often been lost and removal in one piece is generally not possible. These factors often


make it difficult to resect sequestered nuclear material using an intradiscal retrograde approach. In order to achieve adequate decompression, it is therefore necessary to access the extradiscal ventral epidural space directly under continuous visual control.

The most frequent site of lumbar intervertebral disc herniations is in the lower segments. The diameter of the intervertebral foramen decreases in a cranial to caudal direction. Additional narrowing may result from degenerative changes. Particularly at the lower lumbar levels, these anatomic factors often make it difficult to gain extradiscal access to the ventral epidural space under visual control when using the posterolateral approach. Similarly, lateral placement of the endoscope in order to reach the spinal canal tangentially after access has been created is technically difficult due to the preceding passage through soft tissue and the zygapophyseal joint. These problems make it difficult to achieve adequate decompression on a reliable basis when using the posterolateral approach.


For these reasons a lateral transforaminal approach has been developed in recent years.*

In this technique the skin entry point is determined not by measurement in centimeters, but on an individual anatomic basis under radiographic control. The approach permits tangential access to the spinal canal and consequently the direct visualization of the ventral epidural space with continuous irrigation that is required in order to achieve adequate decompression. Used in combination with newly developed endoscopes with a large working channel and corresponding new instruments, shavers, and burrs, this technique has a broad but clearly defined range of indications.

As a guideline for decompression of the spinal canal, caudal and cranial mobility should extend to the middle and start of the pedicle, respectively. Narrowed foramina are no longer a limitation, since they can be broadened. The pelvis can prevent the required lateral access, therefore in an orthograde lateral radiographic view it should not

The skin entry point for the well-known posterolateral approach is measured in centimeters from the midline

With the posterolateral approach the working area is mostly intradiscal

project beyond the middle of the cranial pedicle. At the highest lumbar levels the laterality of the approach is limited by the thoracic and abdominal organs. Because of the greater size of the intervertebral foramen cranially and the possibility of bone resection, a larger radius of action is achieved here and consequently a less lateral access route can be chosen. In the case of intra and extraforaminal decompression operations there are no restrictions. Here, too, a lateral approach is attempted in order to permit atraumatic passage below the exiting spinal root. The surgical access technique for intra- or extraforaminal disc herniations and foraminal stenosis may differ from the conventional technique in order to avoid damaging exiting nerve roots that are dislocated or not definitely localizable. In such cases the extraforaminal approach is used.

Intradiscal procedures such as those required for fusion or infection often call for the posterolateral approach. The approach is always determined by the target point, account being taken of individual pathology and anatomy. Outside


of the established indications for its use, the transforaminal approach has definite limitations.

** Ruetten et al. (2005) An extreme lateral access for the surgery of lumbar disc herniations inside the spinal canal using the full-endoscopic uniportal transforaminal approach. – Technique and prospective results of 463 patients. Spine 30:2570–2578*


Ruetten et al. (2007) Use of newly developed instruments and endoscopes: full-endoscopic resection of lumbar disc herniations via the interlaminar and lateral transforaminal approach. J Neurosurg Spine 6:521-530

The lateral transforaminal approach provides access to the spinal canal in the lower lumbar segments

In the lower lumbar segments the pelvis can prevent the required lateral transforaminal access

The lateral transforaminal approach shifts the working area to the spinal canal

VERTEBRIS lumbar


The full-endoscopic trans- and extraforaminal technique

1. Positioning of the patient

The patient is positioned prone on a radiolucent table with a pelvic and a thoracic roll. Use of a C-arm is required during the operation.


Prone position with pelvic and thoracic roll


2. Determination of lateral access route


Access is determined on the basis of anatomic landmarks under orthograde lateral and posteroanterior fluoroscopic guidance, account being taken of the pathology. Depending on the level, the possibility of injury to abdominal organs must be excluded.

Determination of maximum ventrality on the basis of individual anatomic landmarks and marking of the entry line on the skin


Determination of the intervertebral disc level under orthograde posteroanterior fluoroscopic guidance and determination of the skin entry point

3. Creation of lateral access



After the skin entry point has been determined and the puncture incision made, a spinal cannula is introduced under fluoroscopic guidance, care being taken not to damage neural structures. Positioning in relation to the spinal canal is determined individually in accordance with the target point. The guidewire is then introduced and the spinal cannula removed.

Introduced spinal cannula

The spinal cannula touches the dorsal annulus at the medial pedicular line at the beginning of the spinal canal

The spinal cannula is advanced in the dorsal annulus towards the spinal canal

VERTEBRIS lumbar

The full-endoscopic trans- and extraforaminal technique

Using rotatory movements, the dilator is passed along the guidewire as far as the foramen. After removal of the guidewire it is – depending on the pathology – inserted into the spinal canal.


A beveled working sleeve is then inserted over the dilator and the dilator is removed. Care must be taken to protect neural structures during all working steps.




The guidewire is positioned and the spinal cannula is removed


The dilator is inserted over the guidewire and is in the final position in the spinal canal or dorsal annular defect

The working sleeve is positioned via the dilator and the dilator is removed. The beveled opening is situated inside the spinal canal dorsal to the annulus

VERTEBRIS lumbar

The full-endoscopic **trans- and extraforaminal** technique

4. Performance of operation

The endoscope is passed through the working sleeve. The operation is performed via the intraendoscopic working channel using alternating sets of instruments under full visual control and with continuous irrigation.

The sealing caps for the optic and working sleeve should be used only for brief periods when bleeding obscures vision, since with long operation times and unnoticed obstruction to backflow of irrigation fluid there is a theoretical risk of volume overload and increased pressure within the spinal canal and in the associated and adjacent structures. Experience has shown that, as with all new techniques, the risk of complications is greatest during the learning period.


The lateral approach makes it possible to work in the spinal canal under full visual control

5. Creation of posterolateral access


Obstruction of lateral access by the pelvis or the risk of causing injury to abdominal or thoracic organs in the cranial lumbar segments can necessitate a more posterior or even a posterolateral approach in intradiscal operations. The skin entry point depends on the individual pathology and anatomy and can be either measured in centimeters from the midline or determined by appropriate positioning of the introduced spinal cannula. The subsequent steps, including insertion of the guidewire, the dilator, the sheath, and finally the optic, are the same as in the procedure described above.

Measurement of the skin entry point in centimeters from the midline

The introduced spinal cannula at the desired target point can determine the site of the puncture incision

The maximum laterality of access can be determined on the basis of a preoperative CT scan so as to prevent injury to organs

Operation with posterolateral transforaminal approach


VERTEBRIS lumbar

The full-endoscopic trans- and extraforaminal technique



6. Creation of extraforaminal access

In intra- and extraforaminal intervertebral disk herniations and in foraminal stenosis, the exiting nerve roots may be at increased risk of injury when the access instruments are passed through the foramen. This may necessitate an extraforaminal approach. The skin entry point can be posterolateral to lateral. Instead of being inserted through the foramen into the spinal canal, the spi-

nal cannula is advanced onto the caudal pedicle of the segment to be operated on. This is the safest area in terms of the exiting nerve roots and its use reduces the risk of access-related injury. Subsequently the guidewire, the dilator, and the sheath are likewise advanced until they make bony contact with the pedicle. The anatomic structures of the caudal foramen are then dissected under direct vision, the exiting nerve root is identified, and the operation is performed without damaging the nerve root.


The caudal pedicle is a safe area in terms of the exiting spinal nerve root

Insertion of the spinal cannula as far as the caudal pedicle

Dissection of the anatomic structures of the caudal foramen and the exiting spinal nerve root


7. Bone resection

Resection of bone may be required in order to increase mobility within the spinal canal or when problems arise during access. This can occur, for example, in cases of degenerative or hereditary foraminal stenosis and in recess stenosis operations. The skin entry point can be posterolateral to lateral. After trans- or extraforaminal access has been obtained the bony structures must be dissected for this purpose. In most cases ventral parts of the ascending facet are resected. When parts of the caudal pedicle are resected it must be remembered that this is a weightbearing structure. Extensive resection can weaken the biomechanical structure and result in pedicle fractures.

A range of burrs and bone punches are available for bone resection

Sometimes damaging of the joint in order to reach the medial edge of the ascending facet can not be avoided

In most cases ventral parts of the ascending facet are resected

VERTEBRIS lumbar

The full-endoscopic **trans- and extraforaminal** technique

8. Biportal access

A biportal approach can be required in certain indications such as spondylodiscitis and insertion of implants and when working with special instruments. Access is normally posterolateral using the standard technique. The endoscope can be inserted either unilaterally or in alternating fashion.

Biportal transforaminal access

VERTEBRIS lumbar

The full-endoscopic interlaminar technique

Full-endoscopic interlaminar access

Direct access to the epidural space under continuous vision is a precondition for the performance of satisfactory operations within the spinal canal. When the full-endoscopic transforaminal technique is used, the lateral approach is often required for this purpose. The bony and neural boundaries of the neuroforamen impose limits to mobility and thus also to the indications for operations of this type. Moreover, in the lower lumbar segments the required lateral access can be blocked by the pelvis. It has been found that these limitations make it technically impossible to operate on some pathologies using the full-endoscopic transforaminal approach.

In order to reduce the incidence of surgically induced traumatization of the structures of the spinal canal, it is expedient to make use of anatomically preformed access routes. In addition to the intervertebral foramen, these include the sacral hiatus and the interlaminar window. For technical reasons, epiduroscopy via the sacral hiatus does not permit resection of large structures. This leaves open the possibility of surgical access via the interlaminar window, a long-established and commonly used technique in lumbar spine surgery that was first described in the early 1920s. Various alternative methods were developed in later years, e.g. posterolateral biopsy of the spine in

the late 1940s and intradiscal decompression by means of chemonucleolysis in the early 1970s. Endoscopic inspection of the intervertebral space after open decompression was described in the early 1980s. Full-endoscopic operations were performed mostly using the transforaminal technique with posterolateral access.

A microsurgical technique performed with the aid of a microscope was developed in the late 1970s and went on

Use of the endoscope based on the joystick principle provides mobility

to become the gold standard for interlaminar decompression of the spinal canal. An endoscopically guided technique, or microendoscopic operation, was described in the late 1990s. This used an endoscope to provide visualization of the exposed surgical site on a monitor.

With the conventional technique, the spinal canal has to be opened in order to gain access to the epidural space. This generally involves not only incision of the ligamentum flavum, but also resection of bone. The basic requirement is to achieve adequate access that provides visualization of the spinal canal and permits work with instruments. Problems can arise as a result of traumatization of the access pathway, resection of stabilizing structures,

and – especially in relation to the possible need for revision operations – scarring. The basic role of the microscope is to reduce the size of the access route and provide optimal conditions of light and vision. Resection of structures of the spinal canal is generally unavoidable. The microendoscopic technique provides access with less trauma than does the microscopic technique. Its advantage lies in the smaller distance between the working area and the visualizing system. Visual conditions and illumination are generally worse. It is not a full-endoscopic technique in the strict sense. Nowadays microendoscopic access is sometimes combined with a microscopic surgical technique. Common to all of these techniques is the fact that the access pathway generally has to be bigger than would be necessary for the actual work to be performed in the spinal canal.

The full-endoscopic interlaminar approach was therefore developed in recent years in order to exploit the known advantages of transforaminal operations and arthroscopy.*

* Ruetten et al. (2006) A new full-endoscopic technique for the interlaminar operation of lumbar disc herniations using 6 mm endoscopes: Prospective 2-year results of 331 patients. *Minim Invasive Neurosurgery* 49:80-87

Ruetten et al. (2007) Use of newly developed instruments and endoscopes: full-endoscopic resection of lumbar disc herniations via the interlaminar and lateral transforaminal approach. *J Neurosurg Spine* 6:521-530


VERTEBRIS lumbar

The full-endoscopic interlaminar technique

The fact that the lighting and imaging system with its 25° angle of vision is situated right in the working area makes it possible to minimize traumatization not only of the access pathway, but also of structures of the spinal canal. Continuous irrigation provides excellent visual conditions. Mobility is achieved by handling of the new endoscope using the joystick technique. Neural structures are protected by use of the beveled operating sheath as a nerve hook. When used in conjunction with the newly developed instruments, this is a genuinely minimally invasive technique.

The main indications are pathologies situated within the spinal canal. Attention must be paid to the width of the interlaminar window, which if too narrow may prevent free passage of the endoscope. If this occurs the bone can be resected with a burr until the target point is reached without opening the ligamentum flavum or damaging the zygapophyseal joints. Bone resection should generally be avoided, though in cases of spinal canal stenosis it may be necessitated by the pathology. The incision in the ligamentum flavum can be limited to a few millimeters, since penetration into the spinal canal is facilitated by the elasticity of the ligament. Mobility to the contralateral side is similar to that with conventional operations. In order to minimize resection of structures of the spinal canal, craniocaudal access via neighboring lumbar segments may be considered. The full-endoscopic interlaminar technique permits selective surgery on pathologies situated within the spinal canal with minimal access-induced traumatization. The transforaminal technique is generally more suitable for intradiscal and intra- or extraforaminal work. The transforaminal approach is subject to

more limitations than the interlaminar approach, but causes less tissue damage. Due to anatomic and pathologic factors, the ratio of transforaminal to interlaminar procedures in clinical practice is about 40 to 60.


The interlaminar approach provides outstanding visualization of the structures of the spinal canal

1. Positioning

The patient is positioned prone on a radiolucent table with a pelvic and a thoracic roll. Use of a C-arm is required during the operation.


Prone position with pelvic and thoracic roll

2. Determination of access route


Access is determined on the basis of anatomic landmarks under posteroanterior fluoroscopic guidance, account being taken of the pathology. The skin incision should be made as far medially in the interlaminar window as possible in order to permit insertion in a lateral direction below the obliquely oriented zygapophyseal joints.


Marking of skin entry point

The skin entry point should be as medial as possible

Entry below the zygapophyseal joints should be made possible


Puncture incision

VERTEBRIS lumbal


The full-endoscopic interlaminar technique

3. Creation of access

After the skin entry point has been determined and the puncture incision made, the dilator is inserted as far as the ligamentum flavum under postero-anterior fluoroscopic guidance. The subsequent procedure is performed under lateral fluoroscopic guidance. The working sleeve with a beveled opening is advanced towards the ligament via the dilator and the dilator is removed.


Insertion firstly of the dilator and then of the sheath to the ligamentum flavum under fluoroscopic guidance




4. Performance of operation

The endoscope is passed through the working sleeve. The operation is performed via the intraendoscopic working channel using alternating sets of instruments under full visual control and with continuous irrigation. Once the ligamentum flavum has been opened, the spinal canal can be entered. Mobility is achieved by handling the endoscope using the joystick technique. The operating sheath with its beveled opening serves as a second instrument and can be rotated so as to protect the neural structures.

The operating sheath with its beveled opening can be rotated so as to serve as a second instrument


VERTEBRIS lumbal

The full-endoscopic interlaminar technique

The sealing caps for the optic und working sleeve should be used only for brief periods when bleeding obscures vision, since with long operation times and unnoticed obstruction to backflow of irrigation fluid there is a theoretical risk of volume overload and increased pressure within the spinal canal and in the associated and adjacent structures. In order to reduce the risk of neurologic

damage particularly in the cranial segments, prolonged and continuous excessive medial retraction of the neural structures with the working sleeve must be avoided or else retraction must be performed on an intermittent basis. Experience has shown that, as with all new techniques, the risk of complications is greatest during the learning period.

Use of the endoscope based on the joystick principle provides mobility


Opening of the ligamentum flavum

View of L5-S1 axilla

Bone can be resected as necessary using the available instruments and burrs


The interlaminar approach makes it possible to work in the spinal canal under visual control

VERTEBRIS lumbal

The full-endoscopic interlaminar technique

5. Bone resection


Resection of bone may be required in order to increase mobility within the spinal canal or when problems arise during access. This can occur, for example, in cases of sequestered disc herniations or small interlaminar window and in recess stenosis operations. After access has been obtained, the bony structures are dissected. It may be useful to start decompression at the caudal end of the descending facet. Depending on the pathology, medial parts of the descending or ascending facet or of the caudal and cranial lamina are then resected.



It may be useful to start decompression at the caudal end of the descending facet

A range of burrs and bone punches that can be introduced through the intraendoscopic working channel are available for bone resection

The extent of bone resection depends on the pathology

Lateral bone resection is performed on the floor of the spinal canal in the working area

VERTEBRIS thoracic

Introduction

Depending on the individual pathology and anatomy, transforaminal and interlaminar procedures can also be performed in the thoracic spine. The principal indication is thoracic intervertebral disc herniations without major spinal cord compression that remains symptomatic despite conservative therapy. Only laterally situated pathologies are generally amenable to operation, since manipulation of the spinal cord must be avoided because of the risk of injury and because lateral transforaminal access is prevented by the thoracic organs. When a transforaminal procedure

the lumbar spine. Compared to the lumbar spine, the thoracic spine is subject to a greater risk of injury to neural and surrounding structures and hence also to technical limitations in terms of access and surgical procedure. In cases that are borderline in terms of anatomy, pathology, or clinical features, a conventional operation may be the only suitable option.



Thoracic disc herniation

is planned, a preoperative CT scan should always be performed in order to determine the exact skin entry point and the possibility of free access to the intervertebral disc. Interlaminar access normally requires resection of bone, since the size of the interlaminar window is generally insufficient, especially lateral to the spinal cord. Operations using either access route can be performed anywhere from the cervicothoracic to the thoracolumbar junction and are performed in the same way as in



VERTEBRIS thoracic

1. The full-endoscopic transforaminal technique

Access is determined with the aid of a preoperative CT scan. Structures to be spared are laterally the lung, medially the spinal cord, and ventrally the vessels. Access may be prevented by anatomic or degenerative bone structures such as ribs, transverse processes, or osteophytes. In general a decidedly posterior approach is required.




In order not to cause injury, the spinal cannula should be inserted parallel to the intervertebral space under postero-anterior fluoroscopic guidance; it should lie strictly caudally in the foramen and on making contact with the disc should be situated exactly between the medial and the lateral pedicular line in the foramen. For added safety the spinal cannula can first be advanced onto the bony structures of the intervertebral joint and then directed ventrally along the bone. After the dilator, the operating sleeve, and the optic have been introduced, particular attention should be paid during the operation to the medially situated spinal cord.

2. The full-endoscopic interlaminar technique



The skin entry point is situated over the intervertebral joint/disc on the medial pedicular line, as in cervical foraminotomy. From this point the spinal canal can be reached without manipulation of the spinal cord.

After the dilator, the operating sleeve, and the optic have been introduced, the size of the interlaminar window is generally found to be insufficient to permit entry into the spinal canal without bone resection. A small amount of burring is therefore performed on the medial side

of the joint facets and if necessary on the cranial and the caudal laminae. The lateral part of the spinal canal must be

accessible as far as the intervertebral disc without need to displace the spinal cord medially. There is no limit to craniocaudal extension.

VERTEBRIS

Basic set, **VERTEBRIS lumbar trans-, extraforaminal**
devised by Dr. Rueffen

Set-Nr. 892101111

Article	Types	pcs.
Panoview Plus discoscope, 25°, WL 207 mm, Ø 6.9 x 5.6 mm, working channel Ø 4.1 mm	89210.1254	1
Conical adapter	8791.751	1
Membrane attachment	8792.451	1
Spinal cannula set, 10 pieces, sterile, WL 150 mm, Ø 1.25 mm	4792.803	1
Dilator, Ø. 6.9 mm	89220.1508	1
Working sleeve with bevel, Ø 8.0 mm, WL 185 mm	89220.1078	1
Extension sleeve, Ø 8.0 mm	89220.1408	1
Working sleeve attachment, Ø 8.0 mm	89200.1008	1
Micro-punch, Ø 2.5 mm, WL 360 mm	8792.671	1
Micro-rongeur with long jaws, Ø 2.5 mm, WL 360 mm	89240.1125	1
Nucleus grasping forceps, Ø 3.0 mm, WL 360 mm	89230.1003	1
Nucleus grasping forceps, Ø 4.0 mm, WL 360 mm	89230.1004	1
Tube shaft punch, Ø 4.0 mm, WL 360 mm	89240.1904	1
Atraumatic dissector, Ø 2.5 mm, WL 350 mm	8792.591	1
Atraumatic dissector, Ø 4.0 mm, WL 350 mm	89250.1004	1
X-Tractor	89230.0000	1
Mallet	8866.956	1
Power control generator	2303.001	1
Power stick M4 handle	8564.121	1
Resector, Ø 4.0 mm, WL 350 mm	89970.1004	1
Oval burr, laterally hooded, Ø 4.0 mm, WL 350 mm	89970.1504	1
Oval burr, frontally hooded, Ø 4.0 mm; WL 350 mm	89970.1514	1
"Surgitron" radiofrequency unit	2343.001/ .002	1
Trigger Flex handpiece, complete	8792.691	1
Trigger Flex bipolar electrodes	4792.6912	1

VERTEBRIS

Basic set, **VERTEBRIS lumbar and thoracic interlaminar**
devised by Dr. Ruetten

Set-Nr. 892102222

Article	Types	pcs.
Panoview Plus discoscope, 25°, WL 165 mm, Ø 6.9 x 5.6 mm, working channel Ø 4.1 mm	89210.3254	1
Conical adapter	8791.751	1
Membrane attachment	8792.451	1
Dilator, Ø 6.9 mm	89220.1508	1
Working sleeve with bevel, Ø 8.0 mm, WL 120 mm	89220.3008	1
Working sleeve attachment, Ø 8.0 mm	89200.1008	1
Micro-punch, Ø 2.5 mm, WL 360 mm	8792.671	1
Micro-rongeur with long jaws, Ø 2.5 mm, WL 360 mm	89240.1125	1
Nucleus grasping forceps, Ø 3.0 mm, WL 360 mm	89230.1003	1
Nucleus grasping forceps, Ø 4.0 mm, WL 360 mm	89230.1004	1
Tube shaft punch, Ø 4.0 mm, WL 360 mm	89240.1904	1
Atraumatic dissector, Ø 2.5 mm, WL 350 mm	8792.591	1
Atraumatic dissector, Ø 4.0 mm, WL 350 mm	89250.1004	1
Power control generator	2303.001	1
Power stick M4 handle	8564.121	1
Resector, Ø 4.0 mm, WL 350 mm	89970.1004	1
Oval burr, laterally hooded, Ø 4.0 mm, WL 350 mm	89970.1504	1
Oval burr, frontally hooded, Ø 4.0 mm, WL 350 mm	89970.1514	1
"Surgitron" radiofrequency unit	2343.001/ .002	1
Trigger Flex handpiece, complete	8792.691	1
Trigger Flex bipolar electrodes	4792.6912	1

VERTEBRIS

Basic set, **VERTEBRIS thoracic transforaminal**
devised by Dr. Ruetten

Set-Nr. 892103333

Article	Types	pcs.
Panoview Plus discoscope, 25°, WL 207 mm, Ø 5.9 x 5.0 mm, working channel Ø 3.1 mm	89210.1253	1
Conical adaptor	8791.751	1
Membrane attachment	8792.451	1
Spinal cannula set, 10 pieces, sterile, WL 150 mm, Ø 1.25 mm	4792.803	1
Dilator, Ø 5.9 mm	8792.764	1
Working sleeve with bevel, Ø 7.0 mm, WL 185 mm	89220.1047	1
Extension sleeve, Ø 7.0 mm	89220.1407	1
Working sleeve attachment, Ø 7.0 mm	89200.1007	1
Micro-punch, Ø 2.5 mm, WL 360 mm	8792.671	1
Micro-rongeur with long jaws, Ø 2.5 mm, WL 360 mm	89240.1125	1
Nucleus grasping forceps, Ø 3.0 mm, WL 360 mm	89230.1003	1
Tube shaft punch, Ø 3.0 mm, WL 360 mm	89240.1903	1
Atraumatic dissector, Ø 2.5 mm, WL 350 mm	8792.591	1
Atraumatic dissector, Ø 3.0 mm, WL 350 mm	89250.1003	1
X-Tractor	89230.0000	1
Mallet	8866.956	1
Power control generator	2303.001	1
Power stick M4 handle	8564.121	1
Resector, Ø 3.0 mm, WL 350 mm	89970.1003	1
Oval burr, laterally hooded, Ø 3.0 mm, WL 350 mm	89970.1503	1
Oval burr, frontally hooded, Ø 3.0 mm, WL 350 mm	89970.1513	1
"Surgitron" radiofrequency unit	2343.001/ .002	1
Trigger Flex handpiece, complete	8792.691	1
Trigger Flex bipolar electrodes	4792.6912	1

VERTEBRIS

Endoscopes, working channel 2.7 mm

Article	Types
 Panoview Plus discoscope, 20°, Ø 5.8 x 5.1 mm, WL 205 mm	8792.411
 Panoview Plus discoscope, 20°, Ø 5.8 x 5.1 mm, MRI-compatible, WL 205 mm	8767.412

Endoscopes, working channel 3.1 mm

Article	Types
 Panoview Plus discoscope, 25°, Ø 5.9 x 5.0 mm, WL 207 mm	89210.1253
 Panoview Plus discoscope, 25°, Ø 5.9 x 5.0 mm, WL 165 mm	89210.3253

Endoscopes, working channel 4.1 mm

Article	Types
 Panoview Plus discoscope, 25°, Ø 6.9 x 5.6 mm, WL 207 mm	89210.1254
 Panoview Plus discoscope, 25°, Ø 6.9 x 5.6 mm, WL 165 mm	89210.3254

Endoscope accessories/attachments

Article	Types
 Sealing cap attachment, incl. 10 rubber caps	8792.452
 Sealing caps Ø up to 2.4 mm, pack of 10	89.00
 Sealing membrane	15 479.006
 Membrane attachment	8792.451
 Tap attachment	8791.951
 Conical adapter	8791.751
 O-rings for sealing between fluid adaptor and endoscope, pack of 10	9500.113
 Plug-on eyepiece funnel for connecting C-mount objectives to endoscope optics with plug-on connection	8885.901
 Drip rejector, pack of 20	89200.1000
 Preparation basket for mechanical preparation and sterilization for discoscopes 89210.xxxx	38044.411
 Preparation basket for mechanical preparation and sterilization for discoscopes 8792.411, 8767.412	38044.111
 Antifogging agent	102.02
 Cleaning brush	6.03

VERTEBRIS



Spinal cannula set

Article	Types
 Spinal cannula set, Ø 1.25 mm, pack of 10, sterile, WL 250 mm	4792.802
 Spinal cannula set, Ø 1.25 mm, pack of 10, sterile, WL 150 mm	4792.803

Dilators

Article	Types
 Dilator, Ø 5.9 mm, 1-channel for working sleeves Ø 7.0 mm	8792.763
 Dilator, Ø 5.9 mm, 2-channel for working sleeves Ø 7.0 mm	8792.764
 Dilator, Ø 6.9 mm, 2-channel for working sleeves Ø 8.0 mm	89220.1508

Working sleeves, Ø 7.0 mm

Article	Types
 Working sleeve with 30° bevel, WL 120 mm	89220.3007
 Working sleeve for foraminoplasty, WL 145 mm	89220.1017
 Working sleeve without window, WL 145 mm	89220.1057
 Working sleeves, basic set, WL 165 mm	89220.1907
 Working sleeve with long elevator lip, WL 165 mm	89220.1117
 Working sleeve with long window, WL 165 mm	89220.1087
 Working sleeve for foraminoplasty, WL 165 mm	89220.1007
 Working sleeve with distally closed window, WL 165 mm	89220.1137
 Working sleeve with dual window, WL 185 mm	89220.1027
 Working sleeve with elevator lip, WL 185 mm	89220.1157
 Working sleeve with long elevator lip, WL 185 mm	89220.1167
 Working sleeve with 30° bevel, WL 185 mm	89220.1047
 Working sleeve with 45° bevel, WL 185 mm	89220.1037
 Working sleeve with bevel, WL 185 mm	89220.1147
 Extension sleeve WL 155 mm	89220.1407

VERTEBRIS

Working sleeves, Ø 8.0 mm

Article	Types
 Working sleeve with 30° bevel, WL 120 mm	89220.3008
 Working sleeve for foraminoplasty, WL 145 mm	89220.1018
 Working sleeve, basic set, WL 165 mm	89220.1908
 Working sleeve with long elevator lip, WL 165 mm	89220.1068
 Working sleeve with dual window, WL 185 mm	89220.1028
 Working sleeve with elevator lip, WL 185 mm	89220.1088
 Working sleeve with long elevator lip, WL 185 mm	89220.1098
 Working sleeve with 30° bevel, WL 185 mm	89220.1078
 Working sleeve with 45° bevel, WL 185 mm	89220.1038
 Extension sleeve, WL 155 mm	89220.1408

Fluid adaptors

Article	Types
 Fluid adaptor, Ø 7.0 mm	89220.1307
 Fluid adaptor, Ø 8.0 mm	89220.1308
 Working sleeve attachment, Ø 7.0 mm	89200.1007
 Working sleeve attachment, Ø 8.0 mm	89200.1008
 Sealing caps, pack of 10	89.03

Trepines

Article	Types
 Trepine, WL 195 mm, Ø 5.9 mm, cutting window Ø 3.0 mm	8792.503
 Trepine, WL 195 mm, Ø 5.9 mm, cutting window Ø 5.3 mm	8792.504
 Trepine, WL 195 mm, Ø 6.9 mm, cutting window Ø 6.3 mm	89260.1108

VERTEBRIS

Accessories

Article	Types
 Small spongiosa funnel, for working sleeve Ø 7.0 mm	89220.1517
 Large spongiosa funnel, for working sleeve Ø 7.0 mm	89220.1527
 Spongiosa ram, for working sleeve Ø 7.0 mm	89220.1507

Shaver system

Article	Types
 Power control	2303.001
 Power stick M4	8564.121

Burrs / rotation knives

Article	Types
 Oval burr, Ø 2.5 mm, WL 350 mm, laterally hooded	8792.312
 Resector, Ø 2.5 mm, WL 350 mm, laterally hooded	8792.313
 Nucleus resector, Ø 3.0 mm, WL 350 mm, laterally hooded	89970.1003
 Oval burr, Ø 3.0 mm, WL 350 mm, laterally hooded	89970.1503
 Oval burr, Ø 3.0 mm, WL 350 mm, laterally and frontally hooded	89970.1513
 Nucleus resector, Ø 4.0 mm, WL 350 mm, laterally hooded	89970.1004
 Oval burr, Ø 4.0 mm, WL 350 mm, laterally hooded	89970.1504
 Oval burr, Ø 4.0 mm, WL 350 mm, laterally and frontally hooded	89970.1514
 Oval burr, Ø 4.5 mm, WL 220 mm, laterally hooded	8792.323
 Oval burr, Ø 4.5 mm, WL 220 mm, laterally hooded	8792.321

Irrigation pump system

Article	Types
 Fluid control	2203.001
 Tube system, disposable, with puncture needle, pack of 10	4170.223
 Tube system, disposable, with Safe Lock, pack of 10	4170.224

VERTEBRIS

HF / RF generators

Article	Types
 Surgitron radiofrequency unit, 4 MHz	2343.001/ .002
 Bipolar generator	2352.001/ .002

RF accessories/electrodes

Article	Types
 Trigger Flex handpiece, complete	8792.691
 Trigger Flex bipolar electrodes, pack of 6	4792.6912

HF instruments, bipolar

Article	Types
 Bipolar micro grasping forceps, Ø 2.6 mm, WL 390 mm	89930.1010
 Bipolar inner part, pack of 3	89930.1001
 Tube shaft	89930.1002
 Handle	89930.1000

HF electrodes, bipolar, Ø 2.0 mm, WL 400 mm

Article	Types
 Ring electrode	8765.613
 Button electrode	8765.621
 Stepped ball electrode	8765.612
 Connecting piece	8765.554

HF electrodes, monopolar, WL 400 mm

Article	Types
 Button electrode, Ø 1.6 mm	823.05
 Button electrode, Ø 2.0 mm	823.06
 Button electrode, Ø 2.6 mm	823.08

VERTEBRIS

Forceps / punches, Ø 2.0 mm

Article	Types
 Double-spoon forceps, WL 360 mm	8793.561
 Grasping forceps, WL 360 mm	8793.621
 Micro-punch, WL 360 mm	8793.661

Forceps / punches, Ø 2.5 mm

Article	Types
 Micro-rongeur, WL 290 mm	89240.2025
 Micro-rongeur, double-action, WL 290 mm	89240.2125
 Micro-punch, WL 290 mm	89240.2225
 Micro-bone punch, WL 290 mm	89240.2325
 Micro-rongeur, WL 360 mm	8792.632
 Micro-rongeur, double-action, WL 360 mm	8792.636
 Micro-punch, WL 360 mm	8792.671
 Micro-rongeur, extended jaws, WL 360 mm	89240.1125
 Nucleus grasping forceps, WL 360 mm	89230.1125

Forceps / punches, Ø 3.0 mm

Article	Types+
 Micro-rongeur, WL 290 mm	89240.3003
 Micro-rongeur, double-action, WL 290 mm	89240.3013
 Micro-punch, WL 290 mm	89240.3023
 Tube shaft punch, WL 290 mm	89240.3903
 Micro-rongeur, WL 360 mm	89240.1003
 Micro-rongeur, double-action, WL 360 mm	89240.1013
 Micro-punch, WL 360 mm	89240.1023
 Nucleus grasping forceps, WL 360 mm	89230.1003
 Spreader, WL 360 mm	89230.1803
 Scissors, WL 360 mm	89240.1703
 Tube shaft punch, WL 360 mm	89240.1903

VERTEBRIS

Forceps / punches, Ø 4.0 mm

Article	Types
 Micro-rongeur, WL 290 mm	89240.3004
 Micro-rongeur, double-action, WL 290 mm	89240.3014
 Micro-punch, WL 290 mm	89240.3024
 Tube shaft punch, WL 290 mm	89240.3904
 Micro-rongeur, WL 360 mm	89240.1004
 Micro-rongeur, double-action, WL 360 mm	89240.1014
 Micro-punch, WL 360 mm	89240.1024
 Nucleus grasping forceps, WL 360 mm	89230.1004
 Micro-rongeur, curved, WL 360 mm	89240.1624
 Tube shaft punch, curved, WL 360 mm	89240.1904
 Micro-rongeur, curved, WL 360 mm	89240.1044
 Micro-punch, Ø 2.5 mm, curved, WL 360 mm (fits in 4 mm working channel)	89240.1034

Forceps / punches, Ø 5.2 mm for use through the working sleeve

Article	Types
 Intradiscal grasping forceps, articulating, WL 210 mm	8792.623
 Intradiscal punch, WL 210 mm	8792.663
 Intradiscal rongeur, conical jaw, WL 210 mm	89240.1052

Assorted forceps / punches / scissors for use through the working sleeve


Article	Types
 Punch, Ø 2.7 mm, WL 210 mm	8792.661
 Scissors, Ø 2.7 mm, WL 240 mm	8792.641
 Grasping forceps, Ø 3.4 mm, WL 240 mm	8792.621
 Punch, Ø 3.4 mm, WL 240 mm	8792.662
 Suction punch, Ø 4.5 mm, WL 240 mm	8792.681
 Rongeur, Ø 4.5 x 4.2 mm, WL 210 mm	8791.601
 Rongeur, Ø 4.5 x 4.2 mm, WL 210 mm	8791.691

VERTEBRIS

Hand-held / accessory instruments, sharp-edged

Article	Types
 Annulotome, Ø 2.5 mm, WL 290 mm	89260.2125
 Bone dissector, Ø 2.5 mm, WL 290 mm	89260.2225
 Rasp, Ø 2.5 mm, WL 290 mm	89260.2325
 Trocar, Ø 2.5 mm, WL 290 mm	89260.2425
 Spoon, Ø 2.5 mm, WL 290 mm	89260.2525
 Curette, Ø 2.5 mm, WL 290 mm	89260.2625
 Rasp, Ø 2.5 mm, WL 350 mm	8792.541
 Trocar, Ø 2.5 mm, WL 350 mm	8792.551
 Spoon, Ø 2.5 mm, WL 350 mm	8792.562
 Annulotome, Ø 2.5 mm, WL 350 mm	8792.581
 Curette, Ø 2.5 mm, WL 350 mm	8792.571
 End-cut burr, Ø 3.0 mm, WL 350 mm	89260.1113
 End-cut burr, Ø 4.0 mm, WL 350 mm	89260.1114

Hand-held / accessory instruments, atraumatic

Article	Types
 Elevator, Ø 2.5 mm, WL 290 mm	89250.2025
 Hook probe, Ø 2.5 mm, WL 290 mm	89250.2125
 Probe, Ø 2.5 mm, WL 290 mm	89250.2225
 Dissector, Ø 2.5 mm, WL 350 mm	8792.591
 Dissector, Ø 3.0 mm, WL 350 mm	89250.1003
 Dissector, Ø 4.0 mm, WL 350 mm	89250.1004

VERTEBRIS

Accessories

Article	Types
 Positioning probe	8791.701
 Instrument grasping forceps	8793.856
 "X-Tractor" withdrawal instruments, complete set	89230.0000
 "X-Tractor" clamping device, small	89230.0003
 "X-Tractor" clamping device, large	89230.0004
 "X-Tractor" handle	89230.0008
 Mallet	8866.956

Suction and irrigation accessories

Article	Types
 Flushing valve, foot-operated	89870.0000
 Suction connector	89270.1000
 Sucker, Ø 2.5 mm, WL 290 mm	89270.2025
 Sucker, Ø 4.0 mm, WL 340 mm	89270.1004

Literature

RUETTEN S, KOMP M, MERK H, GODOLIAS G

Use of newly developed instruments and endoscopes: full-endoscopic resection of lumbar disc herniations via the interlaminar and lateral transforaminal approach.

J Neurosurg Spine 2007; 6:521-530

RUETTEN S, KOMP M, GODOLIAS G

Lumbar discectomy with the full-endoscopic interlaminar approach using new developed optical systems and instruments.

WSJ 2006;3:148-156

RUETTEN S, KOMP M, GODOLIAS G

New developed devices for the full-endoscopic lateral transforaminal operation of lumbar disc herniations.

WSJ 2006;3:157-165

RUETTEN S, KOMP M, GODOLIAS G

A new full-endoscopic technique for the interlaminar operation of lumbar disc herniations using 6 mm endoscopes: Prospective 2-year results of 331 patients.

Minim Invasive Neurosurg 2006, 49: 80-87

RUETTEN S, KOMP M, GODOLIAS G

An extreme lateral access for the surgery of lumbar disc herniations inside the spinal canal using the full-endoscopic uniportal transforaminal approach. – Technique and prospective results of 463 patients.

Spine 2005, 30, 2570-2578

RUETTEN S, KOMP M, GODOLIAS G

Full-endoscopic interlaminar operation of lumbar disc herniations using new endoscopes and instruments.

Orthopaedische Praxis 2005, 10, 527532

RUETTEN S

The full-endoscopic interlaminar approach for lumbar disc herniations. In: Mayer HM (ed) Minimally Invasive Spine Surgery. Springer, Berlin Heidelberg New York, 2005, pp 346355

YEUNG AT

Minimally invasive disc surgery with the Yeung Endoscopic Spine System (YESS). Surg Technol Int 8:267-277, 2000

YEUNG AT

The evolution of percutaneous spinal endoscopy and discectomy: state of the art. Mt Sinai J Med 67:327-332, 2000


YEUNG AT, TSOU PM

Posterolateral endoscopic excision for lumbar disc herniation: surgical technique, outcome and complications in 307 consecutive cases. Spine 27:722-731, 2002


YEUNG AT, YEUNG CA

Advances in endoscopic disc and spine surgery: foraminal approach. Surg Technol Int 11:255-263, 2003

Notes

Lined area for notes, consisting of 20 horizontal lines.

M E D I C A L


V E T E R I N A R Y

I N D U S T R I A L

COMBIDRIVE EN

The universal drive system for microsurgery
and endoscopy

COMBIDRIVE EN

Contents

COMBIDRIVE EN – The universal system for microsurgery	4
Shaver applications – for arthroscopy of small joints	5
Shaver applications – for endoscopic spine surgery	6
High-speed applications – Brushless electronic motors with INTRA coupling	7
High-speed applications – Angled handpiece for endoscopic spine burrs with INTRA coupling	8
High-speed applications – Handpieces and angled handpieces for microsurgery with INTRA coupling	9
High-speed applications – Tungsten carbide burrs for straight and angled handpieces	10
High-speed applications – Diamond burrs for straight and angled handpieces	11
High-speed applications – Further burrs and drills for straight and angled handpieces	12
High-speed applications – Micro jig saw with INTRA coupling	13
High-speed applications – Oscillating microsaw with INTRA coupling	14
High-speed applications – Osseo scalpel with sagittal cutting edge with INTRA coupling	15
High-speed applications – Angled handpieces for surgical wires with INTRA coupling	16
High-speed applications – Dermatomes with INTRA coupling	17
Accessories	18
Cleaning, instrument care	19

COMBIDRIVE EN

The universal system for microsurgery

COMBIDRIVE EN drive system

COMBIDRIVE EN is a universal, multidisciplinary drive system for microsurgical and endoscopic interventions. The micro-processor-controlled control unit allows the use of various high-speed motors and shaver handpieces for arthroscopic interventions (rotational speed can be set between 500 and 40,000 rpm). The torque or power generated in the instrument can be preset and is also controlled to keep it constant.

The operating and tool parameters can be set, preset and programmed using an easy-to-clean touch-screen display. The control unit has an adjustable roller pump to cool and irrigate the operating area. The system can be activated and the parameters set using a waterproof foot pedal.

COMBIDRIVE EN set	20951.0000
comprising:	Types
Control unit	20951.0201
Power cable, 1.8 m	20950.0301
Foot pedal, vario	20951.0401
Electronic motor, standard (incl. 2 connecting cables)	80950.0003
Tube set, sterile, disposable	40900.0002
Clip set, 12-piece	40900.0201
Motor handpiece rest	40900.0401
Spray attachment for handpieces	40900.1201
Spray attachment for motors	40900.1202
Stand for irrigation containers	80900.0000

Technical Data:

Speed range: 500 - 40,000 rpm
 Torque: 1 - 10 Ncm
 Pump: 100 ml/min
 Motor power: 200 W
 Dimensions: 290 x 110 x 260 mm
 Weight: 4.5 kg
 Power supply (switchable): 230 / 115 V a.c.

Shaver applications

for arthroscopy of small joints


In terms of ergonomics and performance, the Micro Stick S1 miniature shaver handpiece is ideal for arthroscopy of small joints.

The miniature handpiece allows speeds of 3000 rpm and attachment of various shaver tools from Ø 2 mm to Ø 4.5 mm.

Micro Stick S1

	Ø for shaver tools	Type
Micro Stick S1	2 - 5.5 mm	89951.0001

Accessories for Micro Stick S1 89951.0001 for arthroscopy of small joints							
		Reusable			Disposable (pack of 5)		
		2 mm Ø	3.5 mm Ø	4.5 mm Ø	2 mm Ø	3.5 mm Ø	4.5 mm Ø
	Resector	89971.0032	89971.0033	89971.0034	49971.0032	49971.0033	49971.0034
	Resector, oval		89971.0113	89971.0114		49971.0113	49971.0114
	Resector, curved convex			89971.0254			49971.0254
	Resector, curved concave			89971.0264			49971.0264
	Ball burr with guard		89971.0503			49971.0503	
	Burr oval, acromionizer					49971.0553	

Shaver applications

for endoscopic spine surgery

The Power Stick M4 universal shaver hand-piece is ideal for endoscopic spine surgery.

Power Stick M4
incl. universal connecting cable

	Ø for shaver tools	Type
Power Stick M4 incl. universal connecting cable	2.5 - 4.5 mm	8564.121

Accessories for Power Stick M4 8564.121 for lumbar spine surgery

		Working length	Ø	Types
	Oval burr, laterally hooded	350 mm	2.5 mm	8792.312
	Resector, laterally hooded	350 mm	2.5 mm	8792.313
	Nucleus resector, laterally hooded	350 mm	3 mm	89970.1003
	Oval burr, laterally hooded	350 mm	3 mm	89970.1503
	Oval burr, laterally and frontally hooded	350 mm	3 mm	89970.1513
	Nucleus resector, laterally hooded	350 mm	4 mm	89970.1004
	Oval burr, laterally hooded	350 mm	4 mm	89970.1504
	Oval burr, laterally and frontally hooded	350 mm	4 mm	89970.1514
	Oval burr, laterally hooded	350 mm	4.5 mm	8792.323
	Oval burr, laterally hooded	350 mm	4.5 mm	8792.321

Accessories for Power Stick M4 8564.121 for cervical spine surgery

		Working length	Ø	Types
	Oval burr, laterally hooded	350 mm	2.5 mm	8792.312
	Oval burr, laterally hooded	350 mm	3 mm	89970.1503
	Oval burr, laterally and frontally hooded	350 mm	3 mm	89970.1513

High-speed applications

Brushless electronic motors with INTRA coupling

The motors used in the **COMBIDRIVE EN** system are state of the art. All the motors provide high performance, are extremely reliable and simple to maintain. This is due, among other things, to the sophisticated brushless design. The motors also have the internationally standardised INTRA interface that allows considerable flexibility in terms of user peripherals. Power is supplied to the motors over special, highly flexible and autoclavable cables.

Electronic motor, standard

intended primarily for non-rotating instruments (saws, dermatomes)

Electronic motor, medium

for longer drill inserts

AUTOCLAVE
134° C / 273° F

	Ø	Types
Electronic motor, standard incl. 2 connecting cables	26 mm	80950.0003
Electronic motor, medium incl. 2 connecting cables	21 mm	80950.0002

	Length	Types
Cable for electronic motor, standard	3 m	80950.0203
Cable for electronic motor, medium		80950.0202

High-speed applications

Angled handpiece for endoscopic spine burrs with INTRA coupling

Angled handpiece

AUTOCLAVE
134° C / 273° F

	Max. revs/min	Gearing	Length	Ø	Type
Endoscopic angles handpiece	40,000 rpm	1:1	110 mm	2.35 mm	82950.1301

Hard metal and diamond burrs for endoscopic spine surgery

Tungsten carbide burrs
Ø 3 mm

Diamond burrs
Ø 3 mm and 4 mm

	Ø	Types
Ball burr pack of 3	3.0 mm	82960.3730
	3.0 mm Diamond	82960.3930
	4.0 mm Diamond	82960.3940

	Types
Outer tube with distal guard	82970.1330
Outer tube open	82970.1340

High-speed applications


Handpieces and angled handpieces for microsurgery with INTRA coupling

To be able to satisfy the demands of the widely differing applications, we have developed an extensive variety of straight and angled handpieces.

The various handpieces differ in shape, gearing, dimensions and according to the intended application.

The **gearing** decides the factor by which the tool slows down as the torque is increased.

The **length** of the handpiece decides the depth that can be reached: In general, it is advisable to keep the handpiece as short as possible, however deeper operating sites or endoscopic applications may require longer handpieces.

Suitable tools are described on pages 10-12

AUTOCLAVE
134° C / 273° F

	Max. revs/min	Gearing	Length	For tools		Types
				Shaft Ø	Length from	
Handpiece, short	40,000 rpm	1:1	90 mm	2.35 mm	44 mm	82950.1001
Handpiece, medium			104 mm		70 - 95 mm	82950.1002
Handpiece, long			127 mm		90 - 125 mm	82950.1003

AUTOCLAVE
134° C / 273° F

	Max. revs/min	Gearing	Length	For tools		Types
				Shaft Ø	Length from	
Angled handpiece	40,000 rpm	1:1	130 mm	2.35 mm	70 - 95 mm	82950.1201
					90 - 125 mm	82950.1202

High-speed applications

Tungsten carbide burrs for straight and angled handpieces

The straight and angled handpieces (page 9) are used to drive a burr, drill or diamond tool etc.

The standardised tool holder (2.35 mm) allows extremely flexible use.

We offer a large selection of tools with various shapes and cutting geometries for extremely varied requirements and applications.

Ball-shaped instruments are the classic shape. They are extremely universal and a large selection is available.


Tungsten carbide burrs

Ø from 0.6 to 7.0 mm

Diamond burrs

Ø from 0.6 to 7.0 mm

AUTOCLAVE
134° C / 273° F

	Ø	Shaft Ø	Total length	Types
Tungsten carbide burr set, comprising: 14 burrs and burr stand	0.6 - 7.0 mm	2.35 mm	70 mm	82960.3400
Tungsten carbide burrs, pack of 10	0.6 mm	2.35 mm	70 mm	82960.3406
	0.8 mm			82960.3408
	1.0 mm			82960.3410
	1.4 mm			82960.3414
	1.8 mm			82960.3418
	2.3 mm			82960.3423
	2.7 mm			82960.3427
	3.1 mm			82960.3431
	3.5 mm			82960.3435
	4.0 mm			82960.3440
	4.5 mm			82960.3445
	5.0 mm			82960.3450
	6.0 mm			82960.3460
	7.0 mm			82960.3470

High-speed applications


Diamond burrs for straight and angled handpieces

	Ø	Shaft Ø	Total length		Types	
			70 mm	95 mm		
Diamond burr set comprising: 15 burrs and drill stand	0.6 - 7.0 mm	2.35 mm	•		82960.3600	
Diamond burr set comprising: 10 burrs and drill stand	1.4 - 7.0 mm			•	82960.3800	
Diamond burr, pack of 10	0.6 mm	2.35 mm	•		82960.3606	
	0.7 mm		•		82960.3607	
	0.8 mm		•		82960.3608	
	1.0 mm		•		82960.3610	
	1.4 mm		•		82960.3614	
					•	82960.3814
	1.8 mm		•		82960.3618	
					•	82960.3818
	2.3 mm		•		82960.3623	
					•	82960.3823
	2.7 mm		•		82960.3627	
					•	82960.3827
	3.1 mm		•		82960.3631	
					•	82960.3831
	3.5 mm		•		82960.3635	
					•	82960.3835
	4.0 mm		•		82960.3640	
					•	82960.3840
4.5 mm	•		82960.3645			
			•	82960.3845		
5.0 mm	•		82960.3650			
			•	82960.3850		
6.0 mm	•		82960.3660			
			•	82960.3860		
7.0 mm	•		82960.3670			
			•	82960.3870		

High-speed applications

Further burrs and drills for straight and angled handpieces

Roller shapes are useful when removing flat surfaces tangentially in an anterior direction as opposed to punctiform operation.

Burr, roller-shaped
Ø from 5.0 to 7.0 mm

AUTOCLAVE
134° C / 273° F

	Ø	Shaft Ø	Total length	Types
Burr, roller-shaped pack of 10	5.0 mm	2.35 mm	70 mm	82960.4450
	6.0 mm			82960.4460
	7.0 mm			82960.4470

The drills by Lindemann are intended for drilling precise cylindrical holes.

Drill by Lindemann
Ø from 1.8 to 2.3 mm


AUTOCLAVE
134° C / 273° F

	Ø	Shaft Ø	Total length	Types
Drill by Lindemann, pack of 10	1.8 mm	2.35 mm	70 mm	82960.4618
	2.1 mm			82960.4621
	2.3 mm			82960.4623

High-speed applications

Micro jig saw with INTRA coupling

Saw attachments for cutting and preparing bony structures.

AUTOGLAVE
134° C / 273° F

	Max. revs/min	Gearing	Type
Micro jig saw with spray nozzle (without saw blades)	15,000 rpm	1:1	82950.1401

	Working length	Types
Saw blade, blade thickness 0.4 mm, pack of 6	10 mm	82960.0210
	17 mm	82960.0217
	25 mm	82960.0225

	Shape	Type
Tungsten carbide rasp attachment	0	82960.0000

For preparation and modelling of bone.



Tungsten carbide rasp attachment

High-speed applications

Micro saw, oscillating with INTRA coupling

For dividing and modelling with a more lateral portal.

AUTOCLAVE
134° C / 273° F


	Max. revs/min	Gearing	Oscillation angle	Type
Micro saw, oscillating with spray nozzle (without saw blades)	15,000 rpm	1:1	25°	82950.1601

	Height	Width	Types
Saw blade, with short shaft, blade thickness 0.4 mm	14 mm	6 mm	82960.1006
	14 mm	10 mm	82960.1010
	14 mm	15 mm	82960.1015

High-speed applications

Osseo scalpel with sagittal cutting edge with INTRA coupling

For dividing and modelling with a more pro-grade portal.

AUTOCLAVE
134° C / 273° F


	Max. revs/min	Gearing	Type
Osseo scalpel with sagittal cutting edge and spray nozzle (without saw blades)	15,000 rpm	1:1	82950.1801

	Height	Width	Types
Saw blade, blade thickness 0.4 mm, pack of 6	10 mm	9 mm	82960.0409
	15 mm	6 mm	82960.0606
	15 mm	10 mm	82960.0610

High-speed applications

Angled handpiece for surgical wires with INTRA coupling

This handpiece is used, in particular, in osteosynthesis of small bones, for example OMF, small joints and similar drilling applications (for example, plates, implants).

AUTOCLAVE
134° C / 273° F

	Max. revs/min	Gearing	Length	Types
Angled handpiece for surgical wires with Ø 0.5 to 1.6 mm (without wires)	2,500 rpm	16:1	105 mm	82950.2001
Spare angled handpiece head				82950.2101

	Length	Ø	Types
Surgical wire, pack of 10	70 mm	0.8 mm	82960.2208
	70 mm	1.0 mm	82960.2210
	150 mm	0.8 mm	82960.2408
	150 mm	1.0 mm	82960.2410
	150 mm	1.2 mm	82960.2412
	150 mm	1.5 mm	82960.2415

High-speed applications

Dermatomes with INTRA coupling

Our dermatomes were developed primarily for plastic reconstruction surgery. They are available with blade widths from 12 to 75 mm.

We recommend that they are used with the standard motor 82950.0003.

AUTOCLAVE
134° C / 273° F

	Blade width	Max. revs/min	Types
Dermatome	12 mm	8,000 rpm	82950.2612
	25 mm		82950.2625
	50 mm		82950.2650
	75 mm		82950.2675

	Width	Types
Spare blades for dermatome, pack of 10	12 mm	42960.1412
	25 mm	42960.1425
	50 mm	42960.1450
	75 mm	42960.1475

	Reduction from 75 mm to	Types
Reducing insert for dermatome base plate	50 mm	82960.1501
	25 mm	82960.1502

Accessories

The autoclavable drill stand provides space for up to 16 tools with normal shafts (Ø 2.35 mm) such as drills and burrs.

AUTOCLAVE
134° C / 273° F

	Type
Drill stand for 16 handpiece instruments	82900.0023

	Types
Spare union nut	82900.0800
Spare spanner	82900.1000
Spray tube attachment	82900.0400

Cleaning and instrument care

	Content	Type
Cleaning spray for care of the handpieces	500 ml	40900.1000

	Shape	Types
Spray attachment	 for handpieces	40900.1201
	 for motors	40900.1202

M E D I C A L

V E T E R I N A R Y

I N D U S T R I A L