



# MEMOKATH™



## MEMOKATH™ 028 PROSTATE

- Stent with expansion in one end. There are two versions; SW version for use with flexible cystoscopes and TW version for use with rigid cystoscopes. Both stent versions expand from CH24 to CH42.

Available in lengths<sup>1</sup> from 30 to 70 mm in intervals of 10 mm. Other lengths available on request.

## MEMOKATH™ 044 URETHRA

- Stent with expansion in one end. This stent is made in one version only; TW version for use with rigid cystoscopes. The stent expands from CH24 to CH42.

Available in lengths<sup>1</sup> from 30 to 70 mm in intervals of 10 mm. Other lengths available on request.

## MEMOKATH™ 045 URETHRA

- Stent with expansion in both ends. There are two versions; TW and TTW, both for use with rigid cystoscopes only. The TW stent expands in both ends from CH24 to CH42. The TTW version expands from CH24 to CH42 in one end and from CH24 to CH33 in the other end.

Available in lengths<sup>2</sup> from 10 to 70 mm in intervals of 10 mm. The TW version is also made in 5 mm length. Other lengths available on request.

## MEMOKATH™ 051 URETER

- Stent with expansion in one end. There are two versions; CW version for retrograde insertion and USDW version for antegrade percutaneous insertion. Both stent versions expand from CH10.5 to CH21.

Available in the lengths<sup>1</sup> 30, 60, 100, 150 and 200 mm. Other lengths available on request.

<sup>1</sup> Total length of stent before expansion

<sup>2</sup> Length of stent between expansions

## HOW TO ORDER

- When ordering a stent:
  - 1) Specify the model (without the "0"), i.e. 28 for a prostate stent.
  - 2) Specify the length, i.e. 60 for a 60 mm stent.
  - 3) Specify the type, i.e. SW for use with a flexible cystoscope. Order code: 2860SW.

Distributed by: